

AA- Allergy, Trees Panel

Test ID: 99996

CPT: 82785 , 86003 x 18

Clinical Significance:

See individual tests

Profile Components:

Total IgE	Willow, t12, IgE	Sycamore, t11, IgE
White Ash, t15, IgE	Sweet Gum ,t211, IgE	Red Cedar, t219, IgE
Elm American, t8, IgE	White Pine, t16, IgE	Birch, t3, IgE
Pecan Tree, t22, IgE	Alder, t2, IgE	
Cedar Mountain,t6, IgE	Hazlenut, t4, IgE	
Walnut, t10, IgE	Beech, t5, IgE	
Maple, Box Elde, t1, IgE	Oak, t7, IgE	
Cottonwood, t14, IgE	Olive, t9, IgE	

Labtech Diagnostics

Container:

Serum Separator Tube (SST®)

Transport Temperature:

Room Temperature

Specimen:

Serum

Specimen Stability:

Room temperature: 28 days
Refrigerated: 28 days
Frozen: 1 year

Reject Criteria:

Days Performed:

Mon-Sat

Collection Instructions: